

CLEANTYPE[®] PRIME SERIES

PANEL EASY / PANEL EASY+ (WITH WALL MOUNT)



SPECIAL FEATURES

- Foldable Control Panel Integrated Medical Keyboard
- Backlight brightness adjustable
- Integrated number pad and touch pad
- Damping mechanism when folding down
- Extraordinary Touch Experience
- Cleaning mode to prevent false triggering during cleaning
- EN60601-1-2 Ed 4 Certification
- One USB Extension Port

SPECIFICATION


Number of keys:	109(US-Layout)	Interfaces:	USB 2.0
Keyswitch technology:	carbon contact	Cable length:	150cm
Switching cycles:	approx.2 mill. (per key)	Operating temperature:	0°C~55°C
Protection level:	IP65	Storage temperature:	-30°C~+70°C
Keyboard material:	silicone	Dimensions:	400(L)x223(W)x28.3(T)mm
Housing material:	aluminum	Certification:	approx.2.5kg
Housing color:	metallic grey	Weight:	EN60601-1-2 Ed 4, CE, FCC, RoHS
Backlight:	blue color,5 level adjustable brightness	Warranty:	2 years

ORDER INFORMATION

Cat. No.	Product description	Mouse tray	Color(keyboard)	Layout
PAL-U11220	CLEANTYPE [®] PRIME OR PANEL EASY+		White	US/DE/JP/KR/THAI
PAL-U21220	CLEANTYPE [®] PRIME OR PANEL EASY		Black	US/DE/JP/KR/THAI
PAL-U11230	CLEANTYPE [®] PRIME OR PANEL EASY+	✓	White	US/DE/JP/KR/THAI
PAL-U21230	CLEANTYPE [®] PRIME OR PANEL EASY	✓	Black	US/DE/JP/KR/THAI

Other layouts, configurations and interfaces on request.
 Customized key function on request (Firmware update via USB)

TECHNICAL DRAWING (IN MM)



EN60601-1-2 Ed 4 certified

EN 60601-1-2 Ed 4 applies to the basic safety and essential performance of Medical Equipment (ME) equipment and ME systems in the presence of electromagnetic disturbances and to electromagnetic disturbances emitted by ME equipment and ME systems.

SHORT OPERATING INSTRUCTIONS



Press and hold lock key for 5s to enable Clean Mode.
Clean Mode is indicated by breathing backlight.



Adjust backlight level with this button.
Five level can be adjusted including off.

CLEANING INSTRUCTION

For further details please refer to below link:
<http://www.gett-group.com/content/service-cleaning-instructions>